



## BARIATRIC PATIENT SEATING

# TILT-RECLINE SEATING



The TRC 750 is the choice for tilt and recline positioning of Bariatric patients. Separate manual hand lever controls enable the reclining back and seat tilt to function individually. The seat surface is a metal pan covered with four-inch fabric covered high density foam and is removable if needed. Front riggings are available in Footrests or an optional Full Width Padded Leg Support and Full Width Footboard. Removable Full Length Armrests are standard. Designed for patients weighing up to 750 pounds (340.2 kg).



Push handle and manual control levers two levers to control back recline and two levers to control seat tilt



Back upholstery features adjustable straps to allow upholstery to conform to user

### Areas of Use:

- Long Term Care
- Respiratory Therapy
- Position Change Therapy
- Patient Room
- Bariatric Care

### Accessories and Replacement Parts

- G27209 High Density Foam Seat, Black Fabric, 3' High
- G27208 High Density Foam Backrest, Black Fabric 3' Depth
- G27220 Full Width Padded Leg Support
- G27225 Full Width Footboard
- G27219 Headrest Positioning Wings
- G10033-1 Swing Away Footrests-Pair



back in upright position  
seat in forward tilt position



**TRC28750**  
SWING-AWAY FOOTREST

reclining back in upright position  
seat in reverse tilt position



reclining back in upright position  
seat in forward tilt position



back in full recline position  
seat in reverse tilt position

**TRC28750-FB**  
FULL WIDTH PADDED LEG SUPPORT & FOOTBOARD

### Specifications:

	overall Width		overall Depth		overall Height		Seat Width		Seat Depth		Seat Height floor to pan		Arm Height		Back Height from seat		Back Angle upright/ full recline	Seat Tilt forward	Seat Tilt reverse	Weight Capacity EVENLY DISTRIBUTED		3" (7.62 cm) HD Foam Seat	3" (7.62 cm) HD Foam Backrest	Swing-Away Footrest	Padded Leg Support	Footboard	UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	degrees	degrees	degrees	lb	Kg	Standard	Optional	Does Not Apply	Active Hyperlink		
TRC28750	49.00	124.46	33.50	85.09	52.00	132.08	28.00	71.12	20.00	50.80	16.25	41.28	28.75	73.03	35.00	88.90	10-36	52	2	750	340.19	Standard	Standard	Standard	Standard	Standard	1/ea
TRC28750-FB	49.00	124.46	33.50	85.09	52.00	132.08	28.00	71.12	20.00	50.80	16.25	41.28	28.75	73.03	35.00	88.90	10-36	52	2	750	340.19	Standard	Standard	Standard	Standard	Standard	1/ea

■ Standard    ■ Optional    — Does Not Apply    🔗 Active Hyperlink  
<sup>1</sup> Odd Spec footnote EG: 3.00 inch (1.02 cm) at base



—GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.  
[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
 Gendron is a registered trademark of GF Health Products, Inc.