

EXAM ROOM SUPPLIES

STAINLESS STEEL  
SPONGE BOWL

**Features**

- 18-8 (304) stainless steel
- Seamless construction
- Autoclavable



	Height		Diameter		Volume		UOM
	inch	cm	inch	cm	oz	ml	
3248	2	5.1	4.5	11.4	12	354.9	1 ea
3249	2.5	6.4	5.75	14.6	32	946.4	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
GrafcO is a registered trademark of GF Health Products, Inc.