



SPECIALITY ROLLATORS/WALKERS

GAITSTER

The Lumex Gaitster Forearm Rollator is designed to offer a solution for walking in a more upright position improving posture, reducing lower back strain and giving you the confidence and independence to keep you active and on the go.

Unlike walkers that force you to hunch over putting painful pressure on your wrists and back, the Gaitster is designed to support you in a more upright position giving you better posture so you have more confidence and less discomfort while walking.



Features

- Unique design that gives user a more upright posture to help reduce body strain and provide optimal control and stability when walking
- Contoured and padded adjustable forearm pads that provide comfort and ideal positioning that are adjustable in length AND angle to give the perfect fit and feel for each user
- Ergonomic handgrips that are easy to grip along with easy to engage brakes at your fingertips to use when stopping your rollator and to engage when sitting
- Handle height easily adjusts 40.2" – 47"
- Flexible, padded back support for comfort
- Wide 18" padded seat with fold-up handle offers added comfort
- Cross brace design allows for side-to-side folding for storage and transport and offers added stability to rollator frame
- Large 9.25" reflective casters in front and 8" casters in back for easy maneuverability and indoor or outdoor use
- Deluxe, zippered removable carry pouch for storage of personal items
- Comes with cane holder on side



Specifications & Options:

		overall Width		overall Depth		Width inside handgrips		minimum Handle Height		maximum Handle Height		Wheel front		Wheel rear	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
LX9000B	Standard	27.4	69.6	30.2	76.71	18.5	47	40.2	102.11	47	119.38	9.25	23.5	8	20.32
LX9000G	Standard	27.4	69.6	30.2	76.71	18.5	47	40.2	102.11	47	119.38	9.25	23.5	8	20.32

		Seat Width		Seat Depth		Seat Height from floor		approximate User Height minimum		approximate User Height maximum		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM	Color	Swatch
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
LX9000B	Standard	18	45.72	11	28	24	60.96	5'4"	162.56	6'2"	187.96	23	10.4	300	136	1/ea	Bondi Blue	
LX9000G	Standard	18	45.72	11	28	24	60.96	5'4"	162.56	6'2"	187.96	23	10.4	300	136	1/ea	Graphite	



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.