



ROLLATORS


SET N' GO HEIGHT ADJUSTABLE ROLLATOR

19" WIDE

The Lumex Set n' Go Wide Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one with an added offering for larger individuals. With its enhanced, wider spacing of 19" in between the handgrips and larger 15.75" seat, the Set n' Go Wide is ideal for those users needing a little more room for comfortable movement and seating. Designed with multiple height settings on the frame, the Set n' Go Wide allows the seat height to be adjusted from a hemi 19.5" to a standard 23.5" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

Features

- Easy height adjustment on the frame allows seat height to be adjusted from 19.5" to 23.5" in 1" increments
- Comfortable, 15.75" wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Ergonomic handgrips offer a comfortable grip
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Removable padded backbar with easy release pull knobs
- 8" wheel for both indoor and outdoor use
- Velcro strap attached secures rollator frame when folded for storage and transport
- Designed to meet the requirements of HCPCS: E0143 and E0156
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



		overall Width		Depth at minimum seat Height		Depth at maximum seat Height		Width inside handgrips		minimum Handle Height		maximum Handle Height		Wheel diameter	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4718B	Blue	26.9	68.33	27.5	69.85	31.5	80.01	19.4	49.28	32	81.28	38	96.52	8	20.32
RJ4718R	Burgundy	26.9	68.33	27.5	69.85	31.5	80.01	19.4	49.28	32	81.28	38	96.52	8	20.32
RJ4718S	Silver	26.9	68.33	27.5	69.85	31.5	80.01	19.4	49.28	32	81.28	38	96.52	8	20.32

		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4718B	Blue	15.75	40.01	12	30.48	19.5	49.53	23.5	59.69	16.5	7.48	350	158.76	1 ea
RJ4718R	Burgundy	15.75	40.01	12	30.48	19.5	49.53	23.5	59.69	16.5	7.48	350	158.76	1 ea
RJ4718S	Silver	15.75	40.01	12	30.48	19.5	49.53	23.5	59.69	16.5	7.48	350	158.76	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
DO NOT use this rollator as a wheelchair or transport device.

graham
field

GF1200125 Rev D22