




ROLLATORS

WALKABOUT STEEL KNOCK DOWN ROLLATOR

The Lumex Walkabout Steel Knock Down Rollator offers a comfortable, traditionally designed, value priced solution to mobility for someone looking for an easy to use walking aid. The frame is made of a highly durable steel with 6" easy to maneuver casters. The Walkabout quickly and easily folds into a compact unit for storage and transport.

Features

- Ergonomic hand grips make them easy to grip
- Adjustable handle height to accommodate broad range of user heights
- Durable construction steel frame
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable, hinged, padded backrest can be folded up and down
- 6" casters for ease of maneuverability
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E01056
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ5500R	Burgundy	15.5	39.37	23.5	59.69	22.5	57.15	14	35.56	12	30.48	20	50.80
RJ5500K	Black	15.5	39.37	23.5	59.69	22.5	57.15	14	35.56	12	30.48	20	50.80

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ5500R	Burgundy	6	15.24	30	76.20	34	86.36	18	8.16	325	147.42	1 ea
RJ5500K	Black	6	15.24	30	76.20	34	86.36	18	8.16	325	147.42	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.

graham
field

GF-1800-08RevC23