



CANES

STANDARD WOODEN CANES

Features

- A classic style made of select hardwood that is satin-smooth to the touch
- Available in a variety of attractive finishes
- Comes complete with rubber tip
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Diameter		Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5180A	Walnut Finish	0.875	2.22	36	91.44	250	113.40	6/cs
5181A	Mahogany Finish	0.875	2.22	36	91.44	250	113.40	6/cs
5182A	Walnut Finish	1	2.54	36	91.44	250	113.40	6/cs
5183A	Mahogany Finish	1	2.54	36	91.44	250	113.40	6/cs
5184A	Walnut Finish, Extra Long	1	2.54	42	106.68	250	113.40	6/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.