



CANES

WOOD GRAIN
ALUMINUM CANE

Features

- Stylish, attractive faux wood finish
- Derby-style handle is contoured for comfort
- Black, soft grip provides cushioning for added comfort to reduce hand stress and fatigue
- Easy push-button height adjustment from 29" – 38" with locking collet for added security
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5190WG	Wood Grain	29	73.66	38	96.52	300	136.1	6/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.