



## CANES

# ADJUSTABLE OFFSET CANES

## NITRILE GRIP

### Features

- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Available in a variety of finishes
- Includes laminated color operating instructions
- Soft, warm and contoured for maximum comfort
- Sculpted grip reduces the chance of slippage, even when wet
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5941A	Bronze	31	78.74	39	99.06	250	113.4	6/cs
6333A	Bronze, Extra Long	34.5	87.63	42.5	107.95	250	113.4	1 ea
5942A	Black	31	78.74	39	99.06	250	113.4	6/cs
5943A	Silver	31	78.74	39	99.06	250	113.4	6/cs
5944A	Multi-Pack (2) Bronze, (2) Black, (2) Silver	31	78.74	39	99.06	250	113.4	6/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2015, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF-15-00070RevD22