



## CANES

# FOLDING CANE

### Features

- Cane folds to an 11 inch length for convenient transport and storage
- Derby-style handle is contoured for comfort
- Soft Grip models provide cushioning for added comfort to reduce hand stress and fatigue
- Easy push-button height adjustment from 32.5" - 36.5" with locking collet for added security
- Attractive black frame finish resists marring and scratching
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty



		minimum Height		maximum Height		folded Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	
5950A-4	Black	32.5	82.55	36.5	92.71	11	27.94	250	113.40	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.