



CANES

FASHION
OFFSET CANE

Features

- Available in a variety of attractive finishes
- Offset handle places the user's weight directly over the base for maximum balance and control
- Easy push-button height adjustment with locking collet for added security
- Height adjustment 31"- 40"
- Soft, warm and contoured Nitrile grip for maximum comfort
- Sculpted grip reduces the chance of slippage, even when wet
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5950BI	Blue Ice	31	78.74	40	101.6	250	113.4	4/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.