



CANES

STYLUS OFFSET CANES

Features

- Available in a variety of attractive finishes
- Manufactured with sturdy, extruded aluminum tubing
- Offset handle places the user's weight directly over the base for maximum balance and control
- Soft, contoured foam grip for maximum comfort
- Easy push-button height adjustment from 30" - 39" with locking collet for added security
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5950SB	Blue	30	76.2	39	99.06	250	113.4	6/cs
5950SBK	Black	30	76.2	39	99.06	250	113.4	6/cs
5950SL	Lavender	30	76.2	39	99.06	250	113.4	6/cs
5950SR	Red	30	76.2	39	99.06	250	113.4	6/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1200107/Rev.D02