



CANES

FOLDING
BLIND CANES

Features

- Made of lightweight aluminum tubing
- Features a textured plastic handle for a better grip and a convenient carrying strap
- Cane shaft is covered with red and white reflector tape for nighttime visibility
- Folds down 3 times to 13" and 16" for easy transport and storage
- Available in two lengths - 41" & 46"
- Three (3) year limited warranty



	extended Length		folded Length		UOM
	inch	cm	inch	cm	
5960	41	104.14	13	33.02	1 ea
5961	46	116.84	16	40.64	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.