




CANES

AMERICAN PRIDE OFFSET CANES

Features

- Show your American & Military pride spirit with these designs
- Offset handle places the user's weight directly over the base for maximum balance and wrist support
- Soft, foam grip is comfortable for the hands
- Easy push-button height adjustment from 30" - 39" with locking collet for added security
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity
250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5970CA	Camouflage	30	76.2	39	99.06	250	113.4	6/cs
5970FL	Flag	30	76.2	39	99.06	250	113.4	6/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1600139Rev/C22