



## CANES

# SILVER COLLECTION LOW PROFILE QUAD CANES

### Features

- Offset handle places the user's weight over the base for maximum balance and control
- Low profile, 4-point base further enhances stability
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 29" - 38" in 1" increments
- Tool-free left to right hand adjustment is made easily
- Choice of silver, black, or bronze finish
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0105
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		base Depth		base Width		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6140A	Aluminum, Large Base	12	30.48	8	20.32	29	73.66	38	96.52	300	136.08	4/cs
6140BZ	Bronze, Large Base	12	30.48	8	20.32	29	73.66	38	96.52	300	136.08	4/cs
6140BK	Black, Large Base	12	30.48	8	20.32	29	73.66	38	96.52	300	136.08	4/cs
6141A	Aluminum, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	4/cs
6141BZ	Bronze, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	4/cs
6141BK	Black, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012 GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF1200099Rev102z