



## CANES


# SILVER COLLECTION LOW PROFILE QUAD CANES WITH ORTHO-EASE® GRIP

### Features

- Offset handle places the user's weight over the base for maximum balance and control
- Low profile, 4-point base further enhances stability
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 29" - 38" in 1" increments
- Tool-free left to right hand adjustment is made easily
- Choice of silver or bronze finish
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0105
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty

### Ortho-Ease® Grip

- Grip angle is designed to adjust to comfortably fit the user's hand and wrist angle and to minimize stress on the joints
- Ortho-Ease® grip is ideal for arthritics and for those needing a wider grip for greater support
- Comes complete with convenient wrist strap



		base Depth		base Width		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6120A	Aluminum, Large Base	12	30.48	8	20.32	29	73.66	38	96.52	300	136.08	4/cs
6121A	Aluminum, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	4/cs
6121A-1	Aluminum, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	1 ea
6125A	Bronze, Small Base	8	20.32	6	15.24	29	73.66	38	96.52	300	136.08	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2016 GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF-16010024RevC22