



CANES

BARIATRIC QUAD CANE

Features

- Offset handle places the user's weight over the base for balance and control
- Low profile, 4-point large base is designed for maximum stability and walking assistance
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive silver vein finish resists scratching
- Simple height adjustment from 29" – 38" in 1" increments with locking mechanism for added security
- Tool-free left to right-hand adjustment is made easily
- Intended for larger users
- Comfortable, contoured PVC handgrip
- Designed to meet the requirements of HCPCS code: E0105
- Maximum Weight Capacity 500 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty
- Large Base (12" x 8")

Replacement Parts		UOM
6130AT	Replacement Tips	1 pr
6130AK	Replacement Knob	1 pr



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
6130A	Silver Vein	29	73.66	38	96.52	500	226.80	2/cs
6130A-1	Silver Vein	29	73.66	38	96.52	500	226.80	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field