



CANES

# ADJUSTABLE ALUMINUM CANES

**Features**

- Contoured handle is designed to provide a comfortable, secure hand hold
- Attractive aluminum finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 30" to 39" in 1" increments
- Available with a choice of grips, either durable vinyl or soft foam, and a choice of colors shown below
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
6220A	Vinyl Grip, Aluminum	30	76.20	39	99.06	250	113.4	6/cs
6221A	Soft Grip, Aluminum	30	76.20	39	99.06	250	113.4	6/cs
6222A	Soft Grip, Bronze	30	76.20	39	99.06	250	113.4	6/cs
6223A	Soft Grip, Black	30	76.20	39	99.06	250	113.4	6/cs
6224A	Soft Grip Combo Pack, 2 Aluminum, 2 Bronze, 2 Black	30	76.20	39	99.06	250	113.4	6/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2015, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.