



CANES

IMPERIAL  
OFFSET CANE

**Features**

- Heavy duty steel construction
- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive black finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 29" to 38" in 1" increments
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
6334A	Black	29	73.66	38	96.52	600	272.2	4/cs
6334A-1	Black	29	73.66	38	96.52	600	272.2	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.