

CRUTCHES

UNIVERSAL ALUMINUM CRUTCHES

Features

- Made of lightweight, anodized aluminum
- Reinforced center tube and footpiece provide added strength
- Push button adjustment makes changing height easy
- Height adjustment in 1" increments
- Handgrip position is readily adjusted with wing nuts
- Underarm and handgrip pads are comfortable and durable
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0114
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty

Replacement Parts

		UOM
36100	Replacement Tips, Gray	1 pr
361-GRIP	Replacement Handgrips, Gray	1 pr



		minimum Height		maximum Height		minimum Armpad to Grip		maximum Armpad to Grip		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
3610LF-8	Adult, Latex-Free	45	114.30	53	134.62	15	38.10	19	48.26	5' 2" – 5' 10"	1.57 – 1.78	300	136.08	8 pr/cs
3611LF-8	Tall, Latex-Free	53	134.62	61.5	156.21	16	40.64	20	50.80	5' 10" – 6' 6"	1.78 – 1.98	300	136.08	8 pr/cs
3612LF-8	Youth, Latex-Free	37	93.98	46	116.84	9	22.86	13	33.02	4' 6" – 5' 2"	1.37 – 1.57	300	136.08	8 pr/cs
3613C-8	Combo Pack – 4 Adult, 2 Tall, 2 Youth													8 pr/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0800159RevF22