



CRUTCHES

BARIATRIC
IMPERIAL STEEL CRUTCHES

Features

- Heavy duty steel construction
- Push button height adjustment makes changing height easy
- Height adjustment in 1" increments
- Handgrip position is easily adjusted with wing nuts
- Underarm and handgrip pads are comfortable and durable
- Designed to meet the requirements of HCPCS code: E0114
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		minimum Armpad to Grip		maximum Armpad to Grip		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
3614A	Adult	46	116.84	54	137.16	15	38.10	19	48.26	5' 0" – 6' 0"	1.52 – 1.83	600	272.16	1 pr
3615A	Tall	54	137.16	62	157.48	16	40.64	20	50.80	5' 5" – 6' 4"	1.65 – 1.93	600	272.16	1 pr



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.