



## CRUTCHES

# DELUXE FOREARM CRUTCHES

### Features

- Leg and forearm sections can be adjusted independently for a custom fit
- Vinyl-coated arm cuffs can be shaped to the user's arm for extra comfort
- Contoured vinyl handgrip is durable and comfortable
- Telescoping parts have internal bushings and external lock nuts to silence any noise
- Extra-large, non-skid crutch tips provide added traction and stability
- Constructed with high strength aluminum tubing
- Easy push-button height adjustment with locking collet for added security
- Grip-to-floor height adjusts in 1" increments
- Arm cuff to grip height adjusts in 1" increments
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0110
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		user Height Range		Weight Capacity		UOM
		inch	cm	inch	cm	feet	m	lb	Kg	
6350A	Adult	28	71.12	37	93.98	5' 0" – 6' 0"	1.52 – 1.83	300	136.08	1 pr
6350T	Tall	33	83.82	42	106.68	5' 5" – 6' 4"	1.65 – 1.93	300	136.08	1 pr
6350Y	Youth	22	55.88	31	78.74	4' 10" – 5' 8"	1.47 – 1.73	300	136.08	1 pr



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.