



FOLDING WALKERS

IMPERIAL COLLECTION DUAL RELEASE FOLDING WALKERS X-WIDE

Features

- Increased width and depth accommodates larger users
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0148
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	5" Fixed Wheels No wheels	Width inside handgrips		Width inside rear legs		overall Width		Depth front to rear		Depth folded		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
604070A	●	20	50.80	23	58.42	26.13	66.36	18	45.72	4	10.16	31	78.74	38	96.52	5' 4" - 6' 2"	1.63 - 1.88	600	272.16	2/cs
604070W	●	20	50.80	23	58.42	26.13	66.36	18	45.72	4	10.16	31	78.74	38	96.52	5' 4" - 6' 2"	1.63 - 1.88	600	272.16	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0800176Rev/D22