



FOLDING WALKERS

SINGLE RELEASE FOLDING WALKERS

STANDARD RUBBERIZED PVC GRIP

Features

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Adjustable height fits a broad range of user heights
- Accommodates most Lumex accessories and leg attachments
- Includes laminated color operating instructions
- Handgrip is contoured and textured to provide a better gripping surface
- Designed to meet the requirements of HCPCS code: E0135
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



Handgrip is contoured and textured to provide a better gripping surface

| | | Width | | Width | | overall | | Depth | | Depth | | minimum | | maximum | | user | | Weight | | UOM |
|---------|-------|-------|-------|-------|-------|---------|-------|-------|-------|-------|-------|---------|-------|---------|--------|---------------|-------------|--------|--------|------|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | feet | m | lb | Kg | |
| 613070A | Adult | 17.25 | 43.82 | 19.75 | 50.17 | 23.25 | 59.06 | 19.5 | 49.53 | 4.38 | 11.13 | 32.5 | 82.55 | 39.5 | 100.33 | 5' 4" - 6' 2" | 1.63 - 1.88 | 300 | 136.08 | 4/cs |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1000086RevD22