



FOLDING WALKERS

DUAL RELEASE FOLDING WALKERS

Features

- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Contoured, textured handgrips offer a safe, secure hand hold
- Dual release folding mechanism enables user to fold walker sides independently
- Adjustable height fits a broad range of user heights
- Accommodates all leg attachments
- Folded depth 4"
- Designed to meet the requirements of HCPCS code: (without wheels) E0135 (with wheels) E0143
- Includes laminated color operating instructions
- Models with wheels come with rear glide caps to facilitate use on all surfaces
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		5" Fixed Wheels		No wheels		Width		Width		overall Width		Depth		Depth		minimum Height		maximum Height		user Height Range		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg			
616270A-4	Adult	•	16.5	41.91	20.125	51.12	26	66.04	16.5	41.91	4	10.16	34.5	87.63	41.5	105.41	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs			
616370A-4	Junior	•	16.5	41.91	18.625	47.31	25.38	64.45	15.5	39.37	4	10.16	29.5	74.93	34.5	87.63	4' 7" - 5' 3"	1.4 - 1.6	300	136.08	4/cs			
616070A	Adult	•	16.5	41.91	20.125	51.12	22.75	57.79	16.5	41.91	4	10.16	32.5	82.55	41.5	105.41	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs			



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-15-00016RevD22