



## FOLDING WALKERS

# UPRISE ONYX FOLDING WALKER

### Features

- Combination of folding walker and rising aid all in one
- Secondary handles provide stable assistance from a seated to a standing position
- Attractive, upscale Onyx finish
- Single release folding mechanism is designed to aid users with limited hand dexterity
- Can be used as a portable toilet safety frame
- Durable aluminum tubing provides strength while remaining lightweight
- Adjustable height fits a broad range of user heights
- Contoured, plastic grips for enhanced comfort and a secure hold
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0135
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	Width inside handgrips		Width inside rear legs		overall Width		Depth front to rear		Depth folded		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg		
700175C-2	18	45.72	23	58.42	26	66.04	20	50.8	4.5	11.43	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	400	181.44	2/cs	
700175CR	Unassembled, Retail	18	45.72	23	58.42	26	66.04	20	50.8	4.5	11.43	32	81.28	39	99.06	5' 4" - 6' 2"	1.63 - 1.88	400	181.44	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF0700101RevE22