

## MOBILITY ACCESSORIES

# REPLACEMENT TIPS FOR QUAD CANES

### Features

- 9026 fits large base quad canes: 6118A and 6120A
- 9029 fits small base quad canes: 6119A and 6121A
- Gray

*These tips are recommended for use on Lumex Canes. (We cannot guarantee the performance of these accessories on other industry products.)*



		Tube diameter		UOM
		inch	cm	
9026	Large Base Quad Cane Replacement Tips	0.625	1.59	2 pr/pk
9029	Small Base Quad Cane Replacement Tips	0.5	1.27	2 pr/pk



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF2100098RevC22