



MOBILITY ACCESSORIES

FIXED WHEELS

Features

- Converts standard walker into wheeled walker
- Provides forward mobility with no lateral movement
- The large 5" wheels are ideal for navigating difficult terrain such as pavement, grass or carpet
- Now available with eight adjustment holes to accommodate a broader range of user heights
- Fixed wheels attach to the front of walkers only and are recommended for use with Lumex glide brakes
- Comes complete with glide covers
- Kits of Glide Brakes and 3" or 5" fixed wheels also available
- Kits of 3" or 5" fixed wheels and AutoStop Wheels or Glide Brakes also available
- Designed to fit 1" standard walker tubing
- Designed to meet the requirements of HCPCS code: E0155
- Limited lifetime warranty

These wheels are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



603650A



603630A



603950A

| | | Wheel diameter | | Weight Capacity EVENLY DISTRIBUTED | | UOM |
|---------|---|----------------|-------|---------------------------------------|--------|------|
| | | inch | cm | lb | Kg | |
| 603630A | 3" Fixed Walker Wheels | 3 | 7.62 | 300 | 136.08 | 1 pr |
| 603650A | 5" Fixed Walker Wheels | 5 | 12.70 | 300 | 136.08 | 1 pr |
| 603750A | 5" Fixed Walker Wheels With Leg Extension | 5 | 12.70 | 300 | 136.08 | 1 pr |
| 603950A | 5" Fixed Bariatric Walker Wheels | 5 | 12.70 | 600 | 272.16 | 1 pr |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.