



MOBILITY ACCESSORIES

SWIVEL WHEELS

Features

- Provide a full 360 degree turning radius, swivel wheels make it easy to maneuver in confined areas
- Swivel wheels attach to the front of walkers only and are recommended for use with Lumex glide brakes
- 3" wheels are ideal for navigating over flat surfaces, especially indoors
- The large 5" wheels are ideal for navigating difficult terrain such as pavement, grass or carpet
- Designed to fit 1" standard walker tubing
- Comes complete with glide covers
- Designed to meet the requirements of HCPCS code: E0155
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty

These wheels are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



603850A



603830A

		Wheel diameter		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	lb	Kg	
603830A	3" Swivel Walker Wheels	3	7.62	300	136.08	1 pr
603850A	5" Swivel Walker Wheels	5	12.70	600	272.16	1 pr



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.