

MOBILITY ACCESSORIES

REPLACEMENT TIPS FOR ALUMINUM WALKERS

Features

- Fits most walker legs made of 1 1/8" diameter aluminum tubing

These tips are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



		Tube diameter		UOM
		inch	cm	
LX9141	Walker Replacement Tips	1.125	2.86	4/bag



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF2100112RevC22