



## NON- FOLDING WALKERS

# FREEDOM WALKER

### Features

- Encourages independence and self esteem by allowing freedom of movement and exercise
- Restraint free opening/closing gate with a safety knob
- Mesh available in 13 colors
- Adjustable safety belt to help prevent injuries
- Solid base seat with white cushion for extra comfort
- Constructed of the highest quality health care grade PVC
- Maximum Weight Capacity  
300 lb - EVENLY DISTRIBUTED

*Royal Blue mesh will be provided for online orders.*

If you desire another color, please order via Customer Service at 1-800-347-5678.



80500	internal Width		external Width		minimum Seat Height		maximum Seat Height		minimum Arm Height		maximum Arm Height		overall Height		overall Length		Walking Space front of seat to gate		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg			
	18	45.72	30.5	77.47	20.25	51.44	24.75	62.87	28.75	73.03	33.25	84.46	40.75	103.51	35	88.9	15	38.10	300	136.08	2/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF0700156Rev02z