

EXAM ROOM SUPPLIES

FLAT TYPE
INSTRUMENT TRAYS

Features

- Raised edge all around the perimeter
- Accommodates supplies, instruments, and tools
- 18-8 (304) stainless steel with seamless construction



	overall Length		overall Width		overall Height		UOM
	inch	cm	inch	cm	inch	cm	
3261	13.625	34.6	9.75	24.8	0.625	1.6	1 ea
3262	15.125	38.4	10.5	26.7	0.625	1.6	1 ea
3264	19	48.3	12.5	31.8	0.625	1.6	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.