

THERAPY & REHABILITATION

CERVICAL TRACTION SET

Features

- Deluxe over-door type
- Capacity: 20 pounds
- Sturdy, chrome-plated steel construction
- One size fits all
- Complete with instructions and individually boxed

Includes

- Machine washable, foam-padded, heavy duty head halter with metal support rings
- Hook and loop closures
- Over-door clamp
- Pulley
- 12" Adjustable spreader bar
- 8 Ft. Traction rope
- Double-sealed vinyl water bag
- "S" hook



	adjustable Spreader Bar		Traction Rope		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	feet	m	lb	Kg	
GF1871	12	30.48	8	2.44	20	9.07	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.