



SUPPORT SURFACES



FOAM MATTRESS EXTENSION

The Lumex Foam Mattress Extensions provides the ability to convert an 80" mattress to better fit into an 84" bed frame.



Features

- 4" Foam mattress extension fills the gap when using an 80" mattress on an 84" bed
- Hook and loop fastener secures extension to mattress to reduce shifting
- Soft and durable zippered Chemsafe cover is bacteria resistant, flame retardant, fluid retardant, anti-static, non-allergenic, stain resistant, and easy to clean
- Meets 16 CFR 1633 requirements
- Latex Free



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Shipping Weight ¹		Zippered Cover	Meets 16CFR1633	Made in USA	Mattress in months warranty	Cover in months
			inch	cm	inch	cm	inch	cm	lb	kg					
4INEXT	1/ea	E0184	4	10.16	35	88.90	6	15.24	3	1.36	■	■	■	24	12
4INEXT42	1/ea	E0184	4	10.16	42	106.68	6	15.24	3	1.36	■	■	■	24	12

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply ● Active Hyperlink
¹ approximate - as size changes so does weight to a minor degree

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Product Hyperlinks Enabled
+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.