



## PATIENT LIFTS AND SLINGS



# MAXICARE PATIENT LIFT

## PLx00 SERIES BARIATRIC PATIENT LIFTING SYSTEM

When transferring a bariatric patient, safety is key for both the patient and caregiver. The Gendron Patient Lift is designed to offer a safer and robust solution to patient transfers.

### Features

- Heavy-gauge steel construction
- Four-Point spreader bar with rotation tilts to enhance comfort and safety
- The hand pendant controls all functions of the lift including extending and retracting the legs and the opening and closing of the base
- Lifting and lowering along with the left and right rotational movement are also hand pendant operated
- Lower base height maximizes stability and fits under beds
- 24V DC motor allows quiet smooth operation
- Easy access emergency stop button
- Emergency Power Off button available if needed in event of a power loss
- Powder-coated finish
- Optional digital scale available
- Available in 3 options to support 700 lb (317.5 kg), 800 lb (362.8 kg) and 1,000 (453.5 kg) Weight Capacities

### Warranty

- Frame and spreader bar five (5) year limited
- Actuator one (1) year limited



### Specifications

Item	UOM	Overall Width		Overall Length		ext. legs open Base Width minimum		ext. legs open Base Width maximum		spreader bar Height minimum		spreader bar Height maximum		Floor Clearance		Product Weight		Weight Capacity EVENLY DISTRIBUTED		Actuator Frame & Bar warranty	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg	lb	kg	year	month
PL700	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	40.81	103.66	81.94	208.13	4.5	11.43	351	159.21	700	317.51	12	60
PL800	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	43.31	110.01	78.56	199.54	4.75	12.07	354	160.57	800	362.87	12	60
PL1000	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	43.31	110.01	78.56	199.54	4.75	12.07	354	160.57	100	453.59	12	60



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.