

LAMPS

B-SERIES BATTERY-OPERATED UV HAND-HELD LAMP

Features

- Anodized aluminum housing
- Aluminum reflector optimizes UV irradiance for brighter fluorescent response
- Suitable for long wave lamp application
- Highest UV irradiances of any units of its type
- Compact and lightweight: requires 5 C-batteries (not included)
- Can be operated for approximately five hours on fresh batteries
- NiCad rechargeable batteries can be used
- Transparent plastic tube guard to ensure against accidental breakage
- Tube: 4-Watt (BLE-22B)
- Typical Peak UV Intensity at six inches: 730 $\mu\text{W}/\text{cm}^2$
- Made in USA, except bulbs

USED TO DIAGNOSE only



Tube	Typical Peak UV Intensity at 6" (15.2 cm)	overall Length		overall Width		overall Height		UOM
		inch	mm	inch	mm	inch	mm	
2218 4-Watt (BLE-22B)	750	8	20.3	3.25	8.3	2.5	6.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field

GF1300023Rev/A22