

LAMPS

REPLACEMENT BULBS

Replacement bulbs/lamps for a variety of popular diagnostic instruments.



906



909



911

| | | UOM |
|-----|---|------|
| 906 | Adult Laryngoscope Bulb 2.5v, .38a 8-32 Unc 2a Thread 20 Hours Life TI-1 1/2 Lamp Nickel Plated, Vacuum | 6/bx |
| 909 | Examination Bulb 15.1v, 2.26a 10-32 Tpi Thread 100 Hours Life T-3 Lamp Nickel Plated, Halogen, Xenon | 6/bx |
| 911 | Vaginal Specula Illuminator 4.6v, .81a 12-28 Nf Thread 30 Hours Life TI-1 1/2 Lamp Halogen, Xenon | 6/bx |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.