




## SUPPORT SURFACES



# CLASSIC FOAM MATTRESSES

## 316 & 319 SERIES

The Lumex 316 and 319 Series Foam Mattress provides an excellent value for budget-conscious facilities. The Lumex mattresses are specifically developed with patient comfort and pressure injury management in mind.



### Features

Provides pressure reduction and redistribution of stage 1 pressure ulcers.

### Cover Features

- Vapour permeable (0% vapor transmission)
- Fluid resistant
- Infection Control is aided by Antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic
- Anti-static and flame retardant
- Concealed full zipper (319 Series only)
- Easy to clean and maintain

### Mattress Construction

- Dual layer foam construction with convoluted top layer
- Convoluted top layer provides comfort and allows air to flow and cool the sleep surface.
- Designed to meet the requirements of HCPCS Code E0184

### Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Product <sup>1</sup> Weight		Weight Capacity EVENLY DISTRIBUTED		Raised Side Bolsters	Defined Perimeter	Soft Heel Zone	Soft Heel Slope	Pressure Reduction	Zippered Cover	Meets 16 CFR 1633	Made in USA	Mattress in months	Cover in months
			inch	cm	inch	cm	inch	cm	lb	kg	lb	kg										
31680-1633	1/ea	E0184	35	88.90	80	203.20	6	15.24	19	8.62	300	136.08	—	—	—	—	■	—	■	■	24	12
31980-1633	1/ea	E0184	35	88.90	80	203.20	6	15.24	19	8.62	300	136.08	—	—	—	—	■	■	■	■	24	12
31684-1633	1/ea	E0184	35	88.90	84	213.36	6	15.24	19	8.62	300	136.08	—	—	—	—	■	—	■	■	24	12
31984-1633	1/ea	E0184	35	88.90	84	213.36	6	15.24	19	8.62	300	136.08	—	—	—	—	■	■	■	■	24	12
31675-1633	1/ea	E0184	35	88.90	75	190.50	6	15.24	19	8.62	300	136.08	—	—	—	—	■	—	■	■	24	12
31975-1633	1/ea	E0184	35	88.90	75	190.50	6	15.24	19	8.62	300	136.08	—	—	—	—	■	■	■	■	24	12

■ Standard

■ Optional

■ Available

■ Not Available

■ Accessory

— Does Not Apply

■ Active Hyperlink

<sup>1</sup> approximate - as size changes so does weight to a minor degree



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Product Hyperlinks Enabled  
+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.