



100# AUTOMATIC OBT BASE



Features

- Supports up to 100 lbs of evenly distributed weight
- Automatic/infinite stops within height travel range
- Modified "H" base allows use with chair or bed
- Low profile design allows overbed table to be used with limited clearance beds
- Tops are customizable in size, finish, and material

Base Options:		Width	Length	Height Range		Weight Capacity	UOM
Item	Height	inch	inch	inch - low	inch - hi	lb	
A919000	Extended	18.2	29.3	29.9	43	100	1 ea

Made To Order Top Options:

Item	Top Description	3DL	Laminate	Width	Length	Raised Edge	Flat Edge
		material		inch	inch	top profile	
A28B	Laminate Top	—	■	15	30	—	■
A28-3DLS	3DL OBT Top	■	—	15	30	—	■
A28-RES	3DL Raised Edge	■	—	15	30	■	—
OS3DL0BTR	3DL Oversized	■	—	18	35	■	—
OS3DL0BTF	3DL Oversized	■	—	18	35	—	■



Made To Order Top Finish Options:



Ready Ship - Overbed Table Tops & Base Sets:

Item	Extended Base	base Width		base Length		Height Range		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	lb	kg	
A919053	Solar Oak Finish	18.2	46.2	29.3	74.4	29.9 - 43	75.9 - 109.2	100	45.36	1 ea
A919020	Cocobala Finish	18.2	46.2	29.3	74.4	29.9 - 43	75.9 - 109.2	100	45.36	1 ea
A919029	Windsor Mahogany Finish	18.2	46.2	29.3	74.4	29.9 - 43	75.9 - 109.2	100	45.36	1 ea
A919004	Wild Cherry Finish	18.2	46.2	29.3	74.4	29.9 - 43	75.9 - 109.2	100	45.36	1 ea



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved. Basic American Medical Products is a trademark of GF Health Products, Inc.