



DAILY LIVING AIDS

REUSABLE BED AND CHAIR PAD


Designed to provide the best absorbency to protect a mattress or chair from urinary incontinence and other spills. These are washable and reusable for long periods of use.

Features

- Washable and quick drying
- Triple layer protection
- Quilted cotton comfort top layer
- Extra absorbant soaker layer in middle
- Protective vinyl barrier on bottom with knitted backing
- Rounded corners with reinforced stitching for durability
- Absorbency is approximately: 1,500 ml when it's laying flat on a smooth surface

Care Instructions:

- Machine wash hot water
- Rinse
- Tumble dry low heat
- No softener or bleach



Note: patient's weight and positioning on top can alter the absorbency amount along with eventual wear and tear which can also have an effect. Also, it's not recommended that this quantity of liquid or even a portion of it is left under a patient, as the patient's comfort and well being will be affected when positioned on a wet surface for an extended period of time.

		overall Width		overall Length		UOM
		inch	cm	inch	cm	
D0095-1724R-1	Chair Pad	17	43.18	24	60.96	1 ea
D0095-4235R-1	Bed Pad for Full	42	106.68	35	88.9	1 ea
D0095-5435R-1	Bed Pad for Queen	54	137.16	35	88.9	1 ea
D0095-7835R-1	Bed Pad for King	78	198.12	35	88.9	1 ea
D0095	Bed Pad for Partial	29	73.66	35	88.9	1 ea
D0095-3	Bed Pad Pack	29	73.66	35	88.9	3 ea/pck



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-10-00019-Rev/02z