



DAILY LIVING AIDS

REACHERS

Lightweight aluminum reaching aids help avoid bending down and climbing on objects. A squeeze of the trigger closes jaws for a tight grip. Plus a magnetized tip helps to retrieve small metallic items.

Features:

- Reaches into tight spaces such as behind furniture
- Grabs objects just out of reach in high spaces
- Picks objects up off the ground
- Helps avoid bending and climbing
- Magnetized tip for retrieving small metallic items
- Attractive retail package



| | | Reach | | Jaw Width | | Product Weight | | Lifting Weight | | Jaw Mechanism | Folding Mechanism | Dressing Hook | 360° Jaw Rotation |
|------|------------------|-------|-------|-----------|-------|----------------|------|----------------|------|-------------------------------|-------------------|---------------|-------------------|
| | | inch | cm | inch | cm | lb | kg | lb | kg | | | | |
| 5689 | Basic Reacher | 32 | 81.28 | 3 | 7.62 | 0.5 | 0.23 | 2 | 0.91 | Saw tooth surface | ■ | — | — |
| 5690 | Standard Reacher | 32 | 81.28 | 3 | 7.62 | 0.5 | 0.23 | 2 | 0.91 | Suction cups | ■ | — | — |
| 5691 | Deluxe Reacher | 30 | 76.2 | 2.75 | 6.985 | 0.5 | 0.23 | 2 | 0.91 | Saw tooth & anti-slip surface | — | ■ | ■ |



Magnetized Tip



Basic and Standard versions fold for easy storage and transport



Deluxe version has 360° rotation



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.