

ORTHOPEDIC

DELUXE
FOAM CERVICAL COLLAR

Features

- Medium density, 1" thick foam collar covered in natural stockinette
- Contoured design features chin cutout
- Hook and loop closure



		overall Length		overall Thickness		UOM
		inch	cm	inch	cm	
8601 S	Small	12 - 14	30.5 - 35.6	1	2.54	1 ea
8601 M	Medium	14 - 16	35.6 - 40.6	1	2.54	1 ea
8601 L	Large	16 - 18	40.6 - 45.7	1	2.54	1 ea
8601 XL	X-Large	18 - 20	45.7 - 50.8	1	2.54	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.