

BATH ACCESSORIES

TUB-GUARD[®]
BATHTUB SAFETY RAILS

Features

- Textured, contoured handle provides a good gripping surface
- Clamp holds firmly to tub wall without marring tub surface
- Lightweight, compact design is excellent for travel
- Color: White
- Fits tubs with walls 2 ½" - 6 ¼" thick
- Not recommended for use with fiber glass tubs
- Tool-free installation
- Available in two sizes: standard (12" height) and tall (16 ½" height)
- Tall model provides two gripping surfaces to facilitate both entry and exit from tub
- Attractive retail packaging
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	lb	Kg	
69623A	Standard	12	30.48	250	113.40	3/cs
69663A	Tall	16.5	41.91	250	113.40	3/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.