



## BATH LIFTS

# SPLASH<sup>™</sup> BATH LIFT

### An ultra-lightweight, budget-friendly bath lift

With all the basic elements needed to safely maintain personal hygiene, the Splash bath lift has been designed to offer good value for the budget conscious. Weighing an incredibly low 14.3 lbs, Splash is not only the lightest battery-powered bath lift available in the US, its robust design also comfortably lifts users of up to 280 lbs

### A simple, practical bathing solution

#### Fixed back design

- Optimizes legroom by allowing Splash to fit close to the back of a tub

#### Low seat travel

- Maximizes water immersion for improved personal hygiene

#### Smooth movement and solid, stable design

- Offers a safe bathing experience

#### Slim remote with large hand controls for easy operation

- Can be easily operated by users with small hands and/or limited dexterity

#### Battery level indicator with Stop Descent<sup>™</sup> feature

- Monitors battery power, preventing descent when low

#### Simple to clean surfaces with limited crevices

#### Trim and lightweight for easy moving and handling

#### Total weight of only 14.3 lbs

- Offers quick and uncomplicated fitting and removal from the bath in one complete unit

#### Ultra-compact design

- Provides ultimate ease of carrying
- Can be stored in even the smallest of spaces

#### Subtle, unimposing, and unobtrusive design

- Blends discreetly into any bathroom environment
- Barely noticeable when stored in its lowest position in the bath

#### Slim and Modern look

- Provides a fresh, non medical appearance in the home



**LUMEX<sup>®</sup>**  
**splash<sup>™</sup>**  
 bath lift

Lumex<sup>®</sup> is a registered trademark of GF Health Products, Inc. Splash<sup>™</sup> is a trademark of Mountway Limited

Lumex Splash is powered by Mountway, LTD.



	overall Width		overall Height		overall Depth		Seat Width		Seat Depth		Seat Height Range			Lifting Capacity EVENLY DISTRIBUTED		Net Weight excluding battery		shipping Weight		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg		
5033A-1	26.8	68.07	22.8	57.91	20.1	51.05	13	33.02	13.4	34.04	2.8 - 18.1	7.11 - 45.97	280	127.01	14.3	6.49	19.8	8.98	1 ea	



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
 Lumex is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF1100144/RevB22