




## BATH SEATS

# BATH STOOL

### Features

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Tool-free assembly
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
2060A		13.5	34.29	13.5	34.29	13	33.02	13	33.02	13.5	34.29	17.5	44.45	300	136.08	2/cs
2060R	in Retail Package	13.5	34.29	13.5	34.29	13	33.02	13	33.02	13.5	34.29	17.5	44.45	300	136.08	2/cs
2060R-1	in Retail Package	13.5	34.29	13.5	34.29	13	33.02	13	33.02	13.5	34.29	17.5	44.45	300	136.08	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF1000014RevC22