

BATH SEATS

PLATINUM COLLECTION BATH SEATS

**Features**

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Available with or without backrest
- Comes assembled with tool-free back attachment
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7921A-4	with backrest	20	50.80	18	45.72	20	50.80	12	30.48	14	35.56	20	50.80	350	158.76	4/cs
7931A	without backrest	20	50.80	18	45.72	20	50.80	12	30.48	14	35.56	20	50.80	350	158.76	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.