



BATH SEATS

PLATINUM COLLECTION BATH SEATS

RETAIL

Features

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Available with or without backrest
- Available with or without support arms
- Tool Free design
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7921R	with backrest, standard grey	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	3/cs
7921R-1	with backrest, standard grey	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	1 ea
7921RB	with backrest, steel blue	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	3/cs
7921RB-1	with backrest, steel blue	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	1 ea
7931R	without backrest, standard grey	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	4/cs
7931RB	without backrest, steel blue	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	4/cs
7931RB-1	without backrest, steel blue	20	50.80	19	48.26	20	50.80	12	30.48	14	35.56	22	55.88	350	158.76	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0900066RevF22