



BATH SEATS

BARIATRIC BATH SEAT

Features

- Durable blow-molded plastic seat with optional backrest
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Comes assembled with tool-free back attachment
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- Three (3) year limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Back Support height above seat		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7932A	without Backrest	18	45.72	16.75	42.55	17.75	45.09	12	30.48	14	35.56	17	43.18	N/A		600	272.16	2/cs
7932A-1	without Backrest	18	45.72	16.75	42.55	17.75	45.09	12	30.48	14	35.56	17	43.18	N/A		600	272.16	1 ea
7939A	with Backrest	18	45.72	20.5	52.07	17	43.18	16.25	41.28	16	40.64	20	50.80	13	33.02	600	272.16	2/cs
7939A-1	with Backrest	18	45.72	20.5	52.07	17	43.18	16.25	41.28	16	40.64	20	50.80	13	33.02	600	272.16	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.