



BATH SEATS

PADDED KNOCK DOWN BATH SEATS

Features

- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Knock-down design makes transporting easy
- Quick and easy tool-free assembly
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Back Support height above seat		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7944KD-1	Padded Bath Seat	21	53.34	21	53.34	16	40.64	16	40.64	18	45.72	22	55.88	16	40.64	300	136.08	1 ea
7945KD-1	Padded Bath Seat with Arms	23	58.42	21	53.34	16	40.64	16	40.64	18	45.72	22	55.88	16	40.64	300	136.08	1 ea
7946KD-1	Padded Commode Bath Seat	21	53.34	21	53.34	16	40.64	16	40.64	18	45.72	22	55.88	16	40.64	300	136.08	1 ea
7947KD-1	Padded Commode Bath Seat with Arms	23	58.42	21	53.34	16	40.64	16	40.64	18	45.72	22	55.88	16	40.64	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF12000033RevC22