



COMMODOES

3-IN-1 ALUMINUM COMMODOES

Features

- 3-in-1 design is versatile and economical
- Aluminum frame is lightweight and easy to transport
- Sturdy plastic seat snaps-on easily
- Convenient, tool-free assemble
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0163
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
2195A-4		26	66.04	13.5	34.29	14.5	36.83	17.5	44.45	21.5	54.61	18.5	46.99	300	136.08	4 ea/cs
2195KD-2	Knock Down	26	66.04	13.5	34.29	14.5	36.83	17.5	44.45	21.5	54.61	18.5	46.99	300	136.08	2 ea/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700095Rev1B22