



COMMODES

SILVER COLLECTION 3-IN-1 STEEL DROP ARM COMMUNE

Features

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release drop arm mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Removable backbar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame
- All welded steel construction increases strength and stability
- Durable, grey-coated steel frame is attractive and easy to clean
- Cushioned armrests enhance user comfort
- Sturdy plastic seat snaps-on easily
- Seat height readily adjusts from 18 ¾" to 22 ¾", in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0165
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



Replacement Parts

	UOM
9290 Replacement Tips	4/bag

	overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6433A	23.75	60.33	14	35.56	16.75	42.55	18.75	47.63	22.75	57.79	18	45.72	300	136.08	1 ea
6433A-2	23.75	60.33	14	35.56	16.75	42.55	18.75	47.63	22.75	57.79	18	45.72	300	136.08	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-12000-07/RevB22